## TRAM

## 😫 rifton.

Body support: The TRAM features an innovative support system that secures the client with a single buckle. This simplicity enables the caregiver to quickly prepare the client for a safe transfer.

Built-in up/down switch: Raise and lower the TRAM using the single button built into the right steering handle.

Color-coded clips : Caregivers can quickly hook the metal rings onto one of the TRAM's four clips. The clips have been color-coded for identification in protocols for individual clients.

Telescoping lift column: The TRAM can be raised high enough to support a 193 cm walking client and then retract down to 110 cm for compact storage.

Leg expander: Compact and low-profile, the TRAM passes through a 71 cm doorway. Transfer clients in and out of wider wheelchairs by rotating the leg expander.

Low-base TRAM Smaller casters and a low frame profile give the low-base TRAM a clearance of 10.5 cm, making it perfect for rolling into tight spaces such as beneath hospital beds.





## **Optional Accessories**





## 2019 order form

Use dimension chart to select appropriate size.



			K310 Standard (better for gait training)	K320 Low-base (fits under low beds)
* Frame includes one battery, one battery charger, accessories tote and charger power cord (international power cords available on request)	£	Frame with scale	K301	К321 [
	CHOOSE	Frame	К302 🗆	К322 [
* Thigh straps (pair)	CHOOSE	Narrow	K308	К308
		Wide	К338 🗆	К338 [
Switch pole		I	K666	К666 [
Arm supports (pair)		Small arm prompts	K305	К305 [
	£	Large arm prompts	K315	К315 [
	CHOOSE	Arm platforms	K397	К397 [
		New! Arm platforms w/ handgrips	К396 🗆	К396 [
Pelvic support		Small	К306	К306 [
	CHOOSE	Medium	K316	К316 [
	m	Large	К326 🗌	К326 [
Hip positioner	CHOOSE	Small	K667	К667 [
		Large	K668	К668
Hip positioner pad	CHOOSE	Small	K541	К541 [
		Large	K542	К542 [
Additional battery (recommended for institutional use)			K307	К307 [
Additional thigh straps		Narrow single	К318 🗆	К318 [
	윺	Narrow pair	К308 🗆	КЗ08 [
	CHOOSE	Wide single	К328 🗆	К328
		Wide pair	К338 🗆	К338
Directional lock (pair)		1	К319	К324
Removable back belt			К327 🗆	К327 [
Total			K310	K320

\* Required

For all order enquiries please email admin @medix21.co.nz or Tel 04 233 1627

Exclusive to: Medix21 29 Jarden Mile Ngauranga 6035 Wellington

